

YES1

Recombinant Human YES1 Active GST-His

Catalog No.	CRY001	Quantity:	50 µg
Alternate Names:	Yes, c-yes, HsT441, P61-ES		
Description:	Human YES, Amino acids M ₁ -L ₅₄₃ (as in GenBank entry NM_005433)*, N-terminally fused to GST-HIS ₆ -Thrombin cleavage site *Sequence may contain documented polymorphisms Detailed sequence on request		
Concentration:	0.116 µg/µl		
Gene ID:	7525		
Protein Accession No:	NM_005433		
Source:	Baculovirus infected Sf9 cells		
Molecular Weight:	Theoretical MW _{Fusion Protein} : 90,997 Da		
Formulation:	50 mM Tris-HCl, pH 8.0 + 100 mM NaCl + 5 mM DTT + 4 mM reduced glutathione, 20% glycerol		
Purification:	One-step affinity purification using GSH-agarose		
Product Identity:	YES was confirmed as YES by specific Western Blotting using anti YES antibody		
Specific Activity:	29 pmol/µg×min		

Method for determination of K_m value and specific activity:

- Assay conditions:
60 mM HEPES-NaOH, pH 7.5
3 mM MgCl₂
3 mM MnCl₂
3 µM Na-orthovanadate
1.2 mM DTT
2.5 µg / 50 µl PEG_{20,000}
ATP (variable)
Substrate: Poly(Glu,Tyr)_{4:1}
(Sigma P-0275), 1.0 µg / 50 µl
Recombinant YES: 50 ng / 50 µl
- Filter binding assay
MAFC membrane (Millipore)

Storage & Stability: Store in working aliquots at -80°C. **Avoid repeated freeze-thaw cycles.**




Cell Sciences®
480 Neponset Street
Bldg 12A
Canton, MA 02021

Toll Free: 888-769-1246
Phone: 781-828-0610
Fax: 781-828-0542

E-mail: info@cellsciences.com
Website: www.cellsciences.com

Determination of K_m value for ATP:



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

