


Bovine IgG Affinity Purified

Catalog No.	CS208A CS208B	Quantity:	10 mg 50 mg
Alternate Names:	Immunoglobulin gamma		
Description:	Bovine IgG is purified to > 95% by Protein G affinity chromatography. This product is useful for ELISA, Western Blot, and as an isotype control.		
Source:	Bovine serum		
Concentration:	5.0 mg/ml, lot specific		
Specificity:	May be used as an isotype control as it has no known specificity.		
Isotype:	IgG		
Formulation:	0.02 M Sodium Phosphate, 0.1M NaCl, pH 7.4 without preservative		
Purity:	≥ 95% by SDS-PAGE		
Purification:	Protein G affinity chromatography		
Applications:	ELISA, Western blot, Isotype control		
Storage & Stability:	Store at -80°C for at least 1 year. It is recommended to prepare aliquots of the product upon the initial thaw. Avoid repeated freeze-thaw cycles.		

10% Tris-Glycine SDS-PAGE
Non-reduced, Reduced (3ug)



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.