

Serpine1

Glycosylated Mouse PAI-1 (stable mutant)

Catalog No.	CS217A CS217B	Quantity:	0.5 mg 1.0 mg
Alternate Names:	PAI1, PAI-1, Planh1		
Description:	Glycosylated mouse PAI-1 produced in insect cells with a single conservative mutation (Isoleucine 91 to Leucine) to increase stability.		
Concentration:	2.5 mg/ml		
Gene ID:	18787		
Source:	Insect cell culture		
Molecular Weight:	46 kDa		
Formulation:	Frozen Liquid in 0.05M Sodium Acetate; 0.1M NaCl; 1mM EDTA; pH 5.0		
Purity:	>95% by SDS-PAGE analysis		
Specific Activity:	90% active by uPA titration and SDS-PAGE		
Storage & Stability:	When stored at -80°C, product is stable for 3 years from date of delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.