

Mouse Anti-Human Cardiac Troponin T Clone 3B10 mAb

Catalog No.	CS343A CS343B	Quantity:	100 µg 1.0 mg
Description:	Mouse monoclonal antibody to human Cardiac Troponin T (cTnT). IgG fraction purified by immobilized Protein G. Detects cTnT under non-reducing and reducing conditions. Clone 3B10-G11, isotype IgG1 Kappa.		
Source:	Hybridoma cell culture		
Gene ID:	7139		
Purity:	IgG fraction		
Molecular Weight:	160 kDa		
Host:	Mouse		
Immunogen:	Full length native protein (purified) (Human).		
Isotype:	IgG1 Kappa		
Clone:	3B10		
Concentration:	1.0 mg/ml		
Formulation:	Frozen liquid with 0.05 M Sodium Phosphate + 0.1M NaCl + 1 mM EDTA, pH 6.6		
Target:	Recognizes cTnT		
Cross-Reactivity:	Human		
Applications:	ELISA, WB		
Storage & Stability:	When stored at -80°C, product is stable for 5 years from date of delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

