

Human Prolactin ELISA Kit

Strip well format. Reagents for up to 96 tests

Catalog No.	CS352A CS352B	Quantity:	1 x 96 tests 5 x 96 tests
Intended Use:	This human prolactin antigen assay is intended for the quantitative determination of prolactin antigen in human plasma or breast milk.		
Background:	Prolactin (PRL) is a 199 aa, 23kD peptide hormone that is secreted primarily by the pituitary gland in both males and females, though its major roles are in pregnancy and lactation. Prolactin may have a role in breast cancer development, with higher prolactin levels correlating with postmenopausal breast cancer risk.		
Assay Principle:	Human prolactin will bind to the affinity purified capture antibody coated on the microtiter plate. After appropriate washing steps, biotinylated anti-human prolactin primary antibody binds to the captured protein. Excess primary antibody is washed away and bound antibody, which is proportional to the total prolactin present in the samples, is reacted with peroxidase-conjugated streptavidin. Following an additional washing step, TMB substrate is used for color development at 450nm. A standard calibration curve is prepared along with the samples to be measured using dilutions of human prolactin. Color development is proportional to the concentration of prolactin in the samples.		
Reagents Provided:	<ul style="list-style-type: none">◆ 96 well microtiter strip plate 8X12 removable well strips containing affinity purified anti-human prolactin antibody coated on the blocked and dried surface.◆ 10X Wash Buffer: 1 bottle of 50 ml; bring to 1X using DI water◆ Mouse IgG antigen standard: 1 vial of lyophilized standard◆ Peroxidase anti-mouse IgG antibody: 1 vial of lyophilized biotinylated antibody◆ TMB substrate solution: 1 bottle of 10 ml solution		
Storage and Stability:	All kit components must be stored at 4°C. Store unopened plate and any unused microtiter strips in the pouch with desiccant. Reconstituted standards and primary may be stored at -80°C for later use. DO NOT freeze/thaw the standards and primary antibody more than once. All other unused kit components must be stored at 4°C. Kit should be used no later than the expiration date.		



Reagents and Equipment Required:

- 1-channel pipettes covering 0-10 μ l and 200-1000 μ l
- 12-channel pipette covering 30-300 μ l
- Paper towels or kimwipes
- 50 ml tubes, 1.5 ml centrifuge tubes
- 1N H₂SO₄
- DI water
- Plastic containers with lids
- Microtiter plate spectrophotometer operable at 450 nm
- Microtiter plate shaker with uniform horizontally circular movement up to 300 rpm.

Warnings:

Warning – Avoid skin and eye contact when using TMB substrate solution since it may be irritating to eyes, skin, and respiratory system. Wear safety goggles and gloves.

Precautions:

- **DO NOT** mix any reagents or components of this kit with any reagents or components of any other kit. This kit is designed to work properly as provided.
- **DO NOT** pipette reagents by mouth.
- Always pour substrate out of the bottle into a clean test tube. **DO NOT** pipette out of the bottle as you could contaminate the substrate.
- Keep plate covered except when adding reagents, washing, or reading.
- **DO NOT** smoke, drink, or eat in areas where specimens or reagents are being handled.

Preparation of Reagents:

- **TBS:** 0.1M Tris 0.15 M NaCl, pH 7.4
- **Blocking buffer (BB):** 3% BSA in TBS
- **Wash buffer concentrate:** The wash buffer is supplied in a 10X concentrate and must be diluted 1:10 with deionized water for use with the kit.

Specimen Collection:

The assay measures total human prolactin in the 0.1-100 ng/mL range. Samples with human prolactin levels above 100ng/mL should be diluted in blocking buffer before use. Normal plasma should not require dilution before use in this assay. A 1:2 to 1:4 dilution for breast milk is suggested to ensure that resulting values fall within the linear range of the assay.

Assay Procedure:

Perform assay at room temperature. Vigorously shake plate (300rpm) at each step of the assay.

Preparation of Standard:

Reconstitute standard as directed on the vial to give a 500 ng/mL solution.



Dilution table for preparation of human prolactin standard:

Prolactin concentration (ng/mL)	Dilutions
100	800 μ L (BB) + 200 μ L (500ng/mL)
50	500 μ L (BB) + 500 μ L (100ng/mL)
25	500 μ L (BB) + 500 μ L (50ng/mL)
10	600 μ L (BB) + 400 μ L (25ng/mL)
5	500 μ L (BB) + 500 μ L (10ng/mL)
2.5	500 μ L (BB) + 500 μ L (5ng/mL)
1	600 μ L (BB) + 400 μ L (2.5ng/mL)
0.5	500 μ L (BB) + 500 μ L (1ng/mL)
0.25	500 μ L (BB) + 500 μ L (0.5ng/mL)
0.1	600 μ L (BB) + 400 μ L (0.25ng/mL)
0	500 μ L (BB) Zero point to determine background

NOTE: DILUTIONS FOR THE STANDARD CURVE AND ZERO STANDARD MUST BE MADE AND APPLIED TO THE PLATE IMMEDIATELY.

Standard and Unknown Addition:

Remove microtiter plate from bag. Add 100 μ L standards in duplicate and unknowns to wells. Carefully record position of standards and unknowns. Shake plate at 300 rpm for 30 minutes. Wash wells three times with 300 μ L wash buffer. Remove excess wash by gently tapping plate on paper towel or kimwipe.

Primary Antibody Addition:

Add 10 mL of blocking buffer directly to the primary antibody vial and agitate gently to completely dissolve contents. Add 100 μ L to all wells. Shake plate at 300 rpm for 30 minutes. Wash wells three times with 300 μ L wash buffer. Remove excess wash by gently tapping plate on paper towel or kimwipe.

Streptavidin-HRP Addition:

Dilute 1 μ l into 10ml of blocking buffer and mix well. Add 2.5ml of the diluted secondary into 7.5ml of blocking buffer (for a final dilution of 1:40,000) and add 100 μ l to all wells. Shake plate at 300rpm for 30 minutes. Wash wells three times with 300 μ L wash buffer. Remove excess wash by gently tapping plate on paper towel or kimwipe.



Substrate Incubation:

Add 100 μL TMB substrate to all wells and shake plate for 4-10 minutes. Substrate will change from colorless to different strengths of blue. Quench reaction by adding 50 μL of 1N H_2SO_4 stop solution to all wells when samples are visually in the same range as the standards. Add stop solution to wells in the same order as substrate upon which color will change from blue to yellow. Mix thoroughly and read final absorbance values at 450 nm. For best results read plate immediately.

Measurement:



Set the absorbance at 450 nm in a microtiter plate spectrophotometer. Measure the absorbance in all wells at 450 nm. Subtract zero point from all standards and unknowns to determine corrected absorbance (A_{450}).

Assay Calibration:

Plot A_{450} against the amount of human prolactin in the standards. Fit a straight line through the points using a linear fit procedure. The amount of total human prolactin in the unknowns can be determined from this curve.

A typical standard curve.

(EXAMPLE ONLY, DO NOT USE)



Expected Values:

Human plasma varies between the sexes and is considerably higher during pregnancy. Normal values as defined by the NIH are as follows:

Males: 2-18 ng/mL
Nonpregnant females: 2-29 ng/mL
Pregnant women: 10-209 ng/mL

The concentration of prolactin in human breast milk is relatively high in the days immediately after childbirth, peaking at 157 ± 18 ng/mL on the third day following delivery. Prolactin concentrations then fall quickly to 24 ng/mL by the 13th postpartum day.

Disclaimer:

This information is believed to be correct but does not claim to be all-inclusive and shall be used only as a guide. The supplier of this kit shall not be held liable for any damage resulting from handling or from contact with the above product.

**Example of ELISA Kit Plate Layout:
96 Well Plate**

	1	2	3	4	5	6	7	8	9	10	11	12
A	0	0.1ng/mL	0.25ng/mL	0.5ng/mL	1ng/mL	2.5ng/mL	5ng/mL	10ng/mL	25ng/mL	50ng/mL	100ng/mL	
B	0	0.1ng/mL	0.25ng/mL	0.5ng/mL	1ng/mL	2.5ng/mL	5ng/mL	10ng/mL	25ng/mL	50ng/mL	100ng/mL	
C												
D												
E												
F												
G												
H												

Standards: 22 Wells

Samples: 74 Wells

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
Neponset Valley Tech Park
480 Neponset St., Bldg. 12A
Canton, MA 02021

Toll Free: 888 769-1246
Phone: 781 828-0610
Fax: 781 828-0542

E-mail: info@cellsciences.com
Web Site: www.cellsciences.com