


SERPINE1

Recombinant Rat PAI-1 Biotin Active

Catalog No.	CS445A CS445B	Quantity:	500 µg 1 mg
Alternate Names:	PAI, PAI1, PLANH1		
Description:	Recombinant rat PAI-1 is labeled with biotin at primary amines. This wild type fraction typically contains 50 percent active PAI-1.		
Concentration:	1.0 mg/ml		
Gene ID:	24617		
Source:	<i>E. coli</i>		
Molecular Weight:	43.0 kDa		
Formulation:	Liquid in 0.05 M Sodium Acetate + 0.4 M NaCl + 1 mM EDTA, pH 5.0		
Purity:	>95% by SDS-PAGE analysis		
Specific Activity:	This wild type fraction typically contains 50 percent active PAI-1.		
Storage & Stability:	When stored at -80°C product is stable for 3 years from date of delivery. Avoid repeated freeze-thaw cycles.		



1. Rat PAI-1, Biotin Labeled Active Fraction (3 ug) R
2. Rat PAI-1 (3 ug) + Urokinase (15 ug) R
3. Molecular Weight Markers



Avidin-Alkaline Phosphatase Western Blot
1. Rat PAI-1, Biotin Labeled Active Fraction (0.1 ug) R
2. Rat PAI-1 (3 ug) + Urokinase (2 ug) R
3. Molecular Weight Markers

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

