

KCNK2

Potassium Channel Subfamily K member 2

Catalog No.	CSH0608MP CSH0608PR	Quantity:	10 mg 50 µg
Alternate Names:	TREK1, K2P2.1 Potassium Channel, TPKC1		
Description:	Canonical K ⁺ channels are tetrameric and highly selective for K ⁺ , whereas two-pore domain K ⁺ channels form dimers but with similar pore architecture. KCNK2 is a member of the superfamily of potassium channel proteins containing two pore-forming P domains and its expression is particularly high in brain and heart neurons.		
	The receptor is available in the following formats: stable over-expression cell line, membrane preparation, or purified receptor in HEK293 or CHO. Various tagged versions are available.		
Gene ID:	3776		
UniProtKB:	O95069		
Format:	Cell line, membrane preparation, or purified protein		
Source:	HEK 293 or CHO cells		
Characterization:	Expression of chimera was verified by immunostaining. Receptor demonstrates biological activity when tested in a radioligand assay.		
Affinity Tag Options:	Receptor construct: 2 X TwinStrep-Tag at amino-terminus and HIS ₁₀ at carboxy-terminus		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

