

AURKC

Recombinant Human Aurora Kinase C GST-Tagged Active

Catalog No.	CSI10277	Quantity:	10 µg
Alternate Names:	AIE2, AIK3, STK13, AurC, Aurora-C		
Description:	Recombinant human full length Aurora Kinase C with GST tag. Activated <i>in vitro</i> via autophosphorylation.		
Gene ID:	6795		
Source:	Insect cells		
Molecular Weight:	61.7 kDa		
Formulation:	Liquid in 50 mM Tris, pH 7.5 + 150 mM NaCl + 0.5 mM EDTA + 0.02% Triton X-100 + 2 mM DTT + 50% glycerol		
Purity:	≥ 80% as determined by a Coomassie blue-stained SDS-PAGE gel		
Storage & Stability:	Stable for 6 months in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

