

FRK

Recombinant Human Fyn-Related Kinase GST-Tagged

Catalog No.	CSI10285	Quantity:	10 µg
Alternate Names:	GTK, PTK5, RAK		
Description:	Recombinant human full length FRK, GST-tagged. No special measures were taken to activate this kinase.		
Concentration:	0.25 mg/ml		
Gene ID:	2444		
Protein Accession No:	NP_002022		
Source:	Insect cells		
Molecular Weight:	86.3 kDa		
Formulation:	Liquid in 50 mM Tris, pH 7.5 + 150 mM NaCl + 0.5 mM EDTA + 0.02% Triton X-100 + 2 mM DTT + 50% glycerol		
Specific Activity:	182 nmole of phosphate transferred to poly [Glu,Tyr] 4:1 substrate per minute per mg of total protein at 30°C. Activity determined at a final protein concentration of 4 µg/ml.		
Storage & Stability:	Stable for 6 months in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

