

JAK3

Recombinant Human Janus Kinase 3 (aa 781-1124)

GST-Tagged

Catalog No.	CSI10290	Quantity:	10 µg
Alternate Names:	JAKL, L-JAK, LJAK		
Description:	Recombinant human JAK3 catalytic domain (amino acids 781-1124), GST-tagged. No special measures were taken to activate this kinase.		
Gene ID:	3718		
Protein Accession No:	NP_000206		
Source:	Insect cells		
Molecular Weight:	67.4 kDa		
Formulation:	Liquid in 50 mM Tris, pH 7.5 + 150 mM NaCl + 0.5 mM EDTA + 0.02% Triton X-100 + 2 mM DTT + 50% glycerol		
Storage & Stability:	Stable for 6 months in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

