

## PIM2

### Recombinant Human Pim-2 Oncogene GST-Tagged

<b>Catalog No.</b>	CSI11018	<b>Quantity:</b>	10 µg
<b>Gene ID:</b>	11040		
<b>Protein Accession No:</b>	NP_006866		
<b>Description:</b>	Recombinant human full length PIM2, GST-tagged contains additional C-terminal amino acids KGGRADPAFLYKVV. PIM2 is a proto-oncogene which phosphorylates BAD, which partially inhibits its apoptotic function. No special measures were taken to activate this kinase.		
<b>Source:</b>	Insect cells		
<b>Molecular Weight:</b>	63.7 kDa		
<b>Formulation:</b>	Liquid in 50 mM Tris, pH 7.5 + 150 mM NaCl + 0.5 mM EDTA + 0.02% Triton X-100 + 2 mM DTT + 50% glycerol		
<b>Storage &amp; Stability:</b>	Stable for 6 months in working aliquots at -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

