

ITK

Recombinant Human IL2-inducible T-cell Kinase GST-Tagged Active

Catalog No.	CSI12093	Quantity:	10 µg
Alternate Names:	EMT, LYK, PSCTK2		
Description:	Recombinant human full length ITK, GST-tagged. Activated <i>in vitro</i> via auto-phosphorylation.		
Gene ID:	3702		
Protein Accession No:	NP_005537		
Source:	Insect cells		
Molecular Weight:	99.9 kDa		
Formulation:	Liquid in 50 mM Tris, pH 7.5 + 150 mM NaCl + 0.5 mM EDTA + 0.02% Triton X-100 + 2 mM DTT + 50% glycerol		
Storage & Stability:	Stable for 6 months in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

