

## PTK2

### Recombinant Human Protein Tyrosine Kinase 2 GST-Tagged Active

<b>Catalog No.</b>	CSI12105	<b>Quantity:</b>	10 µg
<b>Alternate Names:</b>	Focal Adhesion Kinase 1, FADK, FAK, FAK1		
<b>Description:</b>	Recombinant human full length PTK2 GST-tagged. Activated <i>in vitro</i> by His-tagged SRC.		
<b>Gene ID:</b>	5747		
<b>Protein Accession No:</b>	NP_722560		
<b>Source:</b>	Insect cells		
<b>Molecular Weight:</b>	146.7 kDa		
<b>Formulation:</b>	Liquid in 50 mM Tris, pH 7.5 + 150 mM NaCl + 0.5 mM EDTA + 0.02% Triton X-100 + 2 mM DTT + 50% glycerol		
<b>Storage &amp; Stability:</b>	Stable for 6 months in working aliquots at -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

