

## Recombinant Human CAMK4 with GST Tag

<b>Catalog No.</b>	CSI12157	<b>Quantity:</b>	10 µg
<b>Alternate Names:</b>	Calcium/calmodulin-dependent Protein Kinase IV, Calspermin, CaMK IV, CaMK-GR, CCDPK		
<b>Description:</b>	Recombinant full length human CAMK4 with GST tag		
<b>Gene ID:</b>	814		
<b>Protein Accession No:</b>	NP_001735		
<b>Source:</b>	<i>E. coli</i>		
<b>Molecular Weight:</b>	79.6 kDa		
<b>Formulation:</b>	Liquid in 50 mM Tris, pH 7.5 + 150 mM NaCl + 0.05% mM EDTA + 0.05% Triton X-100 + 2 mM DTT + 50% glycerol		
<b>Storage &amp; Stability:</b>	Stable for 6 months in working aliquots at -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

