

EPHA4

Recombinant Human EPH receptor A4 (aa 616-887)

GST-Tagged

Catalog No.	CSI13399	Quantity:	10 µg
Alternate Names:	HEK8, SEK, TYRO1		
Description:	Recombinant human EPH receptor A4 catalytic domain (aa 616-887) with GST tag. Activated <i>in vitro</i> via auto-phosphorylation.		
Concentration:	0.46 mg/ml		
Gene ID:	2043		
Protein Accession No:	NP_004429		
Source:	Insect cells		
Molecular Weight:	59.7 kDa		
Formulation:	Liquid in 50 mM Tris, pH 7.5 + 150 mM NaCl + 0.5 mM EDTA + 0.02% Triton X-100 + 2 mM DTT + 50% glycerol		
Specific Activity:	230 nmole of phosphate transferred to poly [Glu,Tyr] 4:1 substrate per minute per mg of total protein at 30°C. Activity determined at a final protein concentration of 2 µg/ml.		
Storage & Stability:	Stable for 6 months in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

