

## FER

### Recombinant Human FER (FPS/FES Related) Tyrosine Kinase (aa 541-822) GST-Tagged

<b>Catalog No.</b>	CSI13406	<b>Quantity:</b>	10 µg
<b>Alternate Names:</b>	TYK3		
<b>Description:</b>	Recombinant human FER catalytic domain (amino acids 541-822) GST-tagged. No special measures were taken to activate this kinase.		
<b>Gene ID:</b>	2241		
<b>Protein Accession No:</b>	NP_005237		
<b>Source:</b>	Insect cells		
<b>Molecular Weight:</b>	60.3 kDa		
<b>Formulation:</b>	Liquid in 50 mM Tris, pH 7.5 + 150 mM NaCl + 0.5 mM EDTA + 0.02% Triton X-100 + 2 mM DTT + 50% glycerol		
<b>Purity:</b>	≥ 95% as determined by SDS-PAGE analysis		
<b>Storage &amp; Stability:</b>	Stable for 6 months in working aliquots at -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

