

MARK1

Recombinant Human MAP/Microtubule Affinity-Regulating Kinase 1 GST-Tagged

Catalog No.	CSI13433	Quantity:	10 µg
Alternate Names:	MARK		
Description:	Recombinant human full length MARK1, GST-tagged. MARK1 is a serine/threonine kinase that is involved in cytoskeletal organization and biogenesis. No special measures were taken to activate this kinase.		
Gene ID:	4139		
Protein Accession No:	NP_061120		
Source:	Insect cells		
Molecular Weight:	116.5 kDa		
Formulation:	Liquid in 50 mM Tris, pH 7.5 + 150 mM NaCl + 0.5 mM EDTA + 0.02% Triton X-100 + 2 mM DTT + 50% glycerol		
Storage & Stability:	Stable for 6 months in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

