

NEK11

Recombinant Human Never In Mitosis Gene A-Related Kinase 11 GST-Tagged

Catalog No.	CSI13442	Quantity:	10 µg
Gene ID:	79858		
Protein Accession No:	NP_079076		
Description:	Recombinant human full length NEK11, GST-tagged. NEK11 belongs to the NIMA family of kinases which are involved in DNA replication and genotoxic stress response. No special measures were taken to activate this kinase.		
Source:	Insect cells		
Molecular Weight:	101.6 kDa		
Formulation:	Liquid in 50 mM Tris, pH 7.5 + 150 mM NaCl + 0.5 mM EDTA + 0.02% Triton X-100 + 2 mM DTT + 50% glycerol		
Storage & Stability:	Stable for 6 months in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

