

## PHKG1

### Recombinant Human Phosphorylase Kinase gamma 1 GST-Tagged

<b>Catalog No.</b>	CSI13453	<b>Quantity:</b>	10 µg
<b>Alternate Names:</b>	PHKG		
<b>Gene ID:</b>	5260		
<b>Protein Accession No:</b>	NP_006204		
<b>Description:</b>	Recombinant full length PHKG1, GST-tagged. PHKG1 is a serine/threonine kinase with glycogenolytic regulatory functions. No special measures were taken to activate this kinase.		
<b>Source:</b>	Insect cells		
<b>Molecular Weight:</b>	72.5 kDa		
<b>Formulation:</b>	Liquid in 50 mM Tris, pH 7.5 + 150 mM NaCl + 0.5 mM EDTA + 0.02% Triton X-100 + 2 mM DTT + 50% glycerol		
<b>Storage &amp; Stability:</b>	Stable for 6 months in working aliquots at -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

