

PHKG2

Recombinant Human Phosphorylase Kinase gamma 2 GST-Tagged

Catalog No.	CSI13454	Quantity:	10 µg
Alternate Names:	GSD9C		
Gene ID:	5261		
Protein Accession No:	NP_000285		
Description:	Recombinant human full length PHKG2 with N-terminal GST tag. PHKG2 is a testis/liver isoform of the phosphorylase kinase gamma subunit. No special measures were taken to activate this kinase.		
Source:	Insect cells		
Molecular Weight:	73.9 kDa		
Formulation:	Liquid in 50 mM Tris, pH 7.5 + 150 mM NaCl + 0.5 mM EDTA + 0.02% Triton X-100 + 2 mM DTT + 50% glycerol		
Storage & Stability:	Stable for 6 months in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

