

H3N2 Influenza A Virus Panama 2007/99

| | | | |
|--------------------|-----------|------------------|--------|
| Catalog No. | CSI15843A | Quantity: | 10 µg |
| | CSI15843B | | 50 µg |
| | CSI15843C | | 1.0 mg |

Description: H3N2 is a subtype of the Influenza A Virus. Its name derives from the forms of the two kinds of proteins on the surface of its coat, Hemagglutinin (H) and Neuraminidase (N). H3N2 exchanges genes for internal proteins with other influenza subtypes. H3N2 has tended to dominate in prevalence over H1N1, H1N2, and Influenza B. H3N2 strain descended from H2N2 by antigenic shift, in which genes from multiple subtypes re-assorted to form a new virus. Both the H2N2 and H3N2 strains contained genes from avian influenza viruses.

Allantoic fluid of 10 days old embryonated eggs, inoculated with Influenza A Virus Panama 2007/99

Source: Allantoic Fluid

Formulation: Liquid in STE + 0.1% Sodium Azide + 0.005% Thimerosal.

Purity: >90% according to 7.5% SDS-PAGE.

Immunological Activity: Has been tested with Anti-Influenza A monoclonal antibodies in ELISA. Serological studies of Influenza A Virus immunogen for antibody production.

Storage & Stability: Stable at 4°C for 1 month. Store desiccated at -20°C.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

