

Human Ascites Fluid Liver Disease/Cirrhosis, Single Donor

Catalog No.	CSI19632A	Quantity:	25 ml
	CSI19632C		1 liter

Description: Ascites, also known as peritoneal cavity fluid, peritoneal fluid excess, hydroperitoneum, is an accumulation of fluid in the peritoneal cavity. Although most commonly due to cirrhosis and severe liver disease, its presence can portend other significant medical problems. Diagnosis of the cause is usually with blood tests, an ultrasound scan of the abdomen, and direct removal of the fluid by needle or paracentesis.

Source: Peritoneal fluid from a single donor with liver disease, cirrhosis.

Appearance: Clear to pale yellow.

Formulation: Frozen Liquid - unprocessed, non-sterile, without preservatives

Country of Origin: USA

Storage & Stability: Store at -20°C to -80°C. Recommended upon initial thaw to prepare working aliquots and store at -20°C to -80°C. **Avoid multiple freeze/thaw cycles.**

Infectious Disease Statement: Not tested for infectious diseases. This material should be handled at Bio-Safety Level 2 (BSL 2) as recommended for potentially infectious human serum or blood specimen in the CCD/NIH manual "Biosafety in Microbiological and Biomedical Laboratories", 2009.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com