

PLAU

Mouse Anti-Human UPA ATF region Clone 724-3A1-8 mAb

Catalog No.	CSI19847A CSI19847B	Quantity:	1.0 mg 10 mg
Alternate Names:	ATF, UPA, URK, u-PA, PLAU		
Description:	Mouse Anti-Human UPA ATF region mAb		
Concentration:	10.0 mg/ml		
Specificity:	Recognizes amino terminal fragment region of Urokinase (uPA)		
Host:	Mouse		
Immunogen:	Recombinant full length wild-type human protein (Glycosylated).		
Isotype:	IgG		
Clone:	724-3A1-8		
Formulation:	Liquid in 0.05M Sodium Phosphate + 0.1M NaCl + 1 mM EDTA, pH 6.6		
Purification:	Protein G purification		
Applications:	ELISA Western Blot		
Storage & Stability:	Store at -80°C. Minimum shelf life is 5 years.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.