

Biotin labeled Anti Human 515 kDa LRP mAb

Catalog No.	CSI19867A CSI19867B	Quantity:	100 µg 1.0 mg
Description:	Biotin labeled monoclonal antibody to human 515 kDa heavy chain of LRP. Purified by immobilized Protein A. IgG1k class.		
Purity:	IgG fraction		
Specificity:	Recognizes the 515kDa alpha chain of Low Density LRP/alpha2 MR.		
Host:	Mouse		
Immunogen:	Full length native protein (purified) (Human).		
Isotype:	IgG1 kappa		
Formulation:	Liquid in 0.05M Sodium Phosphate + 0.1M NaCl + 1mM EDTA, pH 6.6		
Purification:	Purified by immobilized Protein A. IgG1k class.		
Cross-Reactivity:	Human		
Applications:	Dot, EM, Flow Cyt, IHC-Fr, IHC-P, IF, WB, BL		
Application Notes:	The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	When stored at -80°C, product is stable for 5 years from date of delivery.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

