

**F9**

## **Rat Anti-Mouse Factor IX & IXa Clone 9041 mAb**

<b>Catalog No.</b>	CSI19869A CSI19869B	<b>Quantity:</b>	100 µg 0.5 mg
<b>Alternate Names:</b>	Cf9, Cf-9, F9		
<b>Description:</b>	Rat Anti-Mouse Factor IX & IXa Clone 9041 mAb		
<b>Concentration:</b>	5.0 mg/ml		
<b>Gene ID:</b>	14071		
<b>Purity:</b>	IgG fraction		
<b>Specificity:</b>	Recognizes Factor IX and IXa		
<b>Host:</b>	Rat		
<b>Immunogen:</b>	Full length, purified native mouse protein.		
<b>Isotype:</b>	IgG		
<b>Clone:</b>	AMIXA-9041		
<b>Formulation:</b>	Liquid in 50% glycerol		
<b>Purification:</b>	IgG fraction		
<b>Applications:</b>	ELISA: Detects Factor IX and Factor IXa. Western Blot: Detects Factor IX and Factor IXa heavy chain reduced and non-reduced Affinity Purification		
<b>Storage &amp; Stability:</b>	Store at -20°C. Stable for 1 year from delivery. <b>Avoid repeated freeze-thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**