

F7

Rat Anti-Mouse Factor VII Clone 9031 mAb

Catalog No.	CSI19870A CSI19870B	Quantity:	100 µg 0.5 mg
Alternate Names:	Cf7, FVII, mfVII, F7		
Description:	Rat Anti-Mouse Factor VII Clone 9031 mAb		
Concentration:	5.0 mg/ml		
Gene ID:	14068		
Purity:	IgG fraction		
Specificity:	Recognizes Factor VII and VIIa		
Host:	Rat		
Immunogen:	Full length, purified native mouse protein.		
Isotype:	IgG		
Clone:	AMVII-9031		
Formulation:	Liquid in 50% glycerol		
Purification:	IgG fraction		
Applications:	ELISA: Detects recombinant Mouse Factor VII and Factor VIIa. Western Blot: Detects recombinant Mouse Factor VII and Factor VIIa non-reduced.		
Storage & Stability:	Store at -20°C. Stable for 1 year from delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.