

F10

Rat Anti-Mouse Factor X Clone 9051 mAb

Catalog No.	CSI19873A CSI19873B	Quantity:	100 µg 0.5 mg
Alternate Names:	fX, Cf10, F10		
Description:	Rat Anti-Mouse Factor X Clone 9051 mAb		
Concentration:	5.0 mg/ml		
Gene ID:	14058		
Purity:	IgG fraction		
Specificity:	Recognizes heavy chain of Factor X		
Host:	Rat		
Immunogen:	Full length, purified native mouse protein.		
Isotype:	IgG		
Clone:	AMX-9051		
Formulation:	Liquid in 50% glycerol		
Purification:	IgG fraction		
Applications:	ELISA: Detects Factor X. Western Blot: Detects Factor X heavy chain non-reduced		
Storage & Stability:	Store at -20°C. Stable for 1 year from delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.