

PIg

Rat Anti-Mouse Plasminogen (Clone 9130) mAb

Catalog No.	CSI19876A CSI19876B	Quantity:	100 µg 0.5 mg
Alternate Names:	Plasminogen, PLG, PLMN		
Description:	Rat Anti-Mouse Plasminogen monoclonal antibody recognizes mouse plasminogen and plasmin.		
UniProt ID:	P20918		
Source:	Rat ascites		
Specificity:	Recognizes mouse plasminogen and plasmin. Does not react with human or bovine plasminogen.		
Immunogen:	Full length, purified native mouse protein		
Purity:	> 95% by SDS-PAGE		
Clone:	9130		
Concentration:	5.0 mg/ml		
Formulation:	Liquid in 50% glycerol		
Purification:	Protein G affinity chromatography		
Applications:	ELISA: Detects plasminogen and plasmin Western Blot: Detects plasmin non-reduced only. Detects plasminogen reduced and non-reduced.		
Storage & Stability:	Store at -20°C to -80°C for up to 1 year. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com