

Chromogenic Substrate for Plasmin

Catalog No.	CSI20043A CSI20043B	Quantity:	25 mg 100 mg
Description:	H-D-Val-Leu-Lys-pNA •2HCl		
Molecular Weight:	551.6 Da		
Formulation:	Lyophilized powder. Each 25 mg of substrate is lyophilized with 60 mg mannitol as a bulking agent.		
Purity:	>95% by SDS-PAGE analysis		
Reconstitution:	Resuspend to 10 mM in 1 mM HCl. Dilute to 0. 2 mM in assay buffer.		
Storage & Stability:	Store at -80°C.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

