

F9

Rat Anti-Mouse Factor IX Clone 942 mAb

Catalog No.	CSI20055A CSI20055B	Quantity:	100 µg 0.5 mg
Description:	Anti-Mouse Factor IX monoclonal antibody detects Factor IX and Factor IXa in ELISA, and Factor IX and Factor IXa non-reduced in western blot.		
Concentration:	5.0 mg/ml		
Gene ID:	14071		
Purity:	IgG fraction		
Specificity:	Mouse Factor IX and IXa		
Host:	Rat		
Immunogen:	Full length purified native mouse protein		
Isotype:	IgG		
Clone:	942		
Formulation:	Rat ascites containing 50% glycerol		
Applications:	ELISA, Western Blot		
Storage & Stability:	Store at -20°C. Stable at 20°C for 1 year from delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

