

## Cynomolgus Monkey IgG Isotype Control pAb

<b>Catalog No.</b>	CSI20163A CSI20163B	<b>Quantity:</b>	10 mg 50 mg
<b>Description:</b>	Cynomolgus Monkey IgG, isotype control, purified using Protein A affinity chromatography.		
<b>Concentration:</b>	≥ 10.0 mg/ml, lot specific		
<b>Specificity:</b>	This antibody has been selected as an isotype control and has no known specificity.		
<b>Host:</b>	Cynomolgus monkey		
<b>Isotype:</b>	IgG		
<b>Formulation:</b>	20 mM Sodium Phosphate, 150 mM NaCl, pH 7.4, preservative free		
<b>Purification:</b>	Purified from normal serum by immobilized Protein A		
<b>Purity:</b>	> 95% by SDS-PAGE and preservative free		
<b>Applications:</b>	ELISA, Western Blot		
<b>Storage &amp; Stability:</b>	Store at -20°C to -80°C for up to 1 year. Upon initial thaw, prepare working aliquots and store at -20°C to -80°C. <b>Avoid repeated freeze/thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**