

FGA, FGB, FGG Native Rabbit Fibrinogen FITC

Catalog No.	CSI20227A CSI20227B	Quantity:	10 mg 40 mg
Description:	Rabbit Fibrinogen labeled at the amino terminus with FITC. Prepared from fresh rabbit plasma using several chromatographic steps. Plasminogen depleted by lysine affinity chromatography.		
UniProt ID:	P14461, P14478		
Gene ID:	100353616, 100351341, 10035159		
Source:	Rabbit plasma		
Molecular Weight:	330 kDa		
Formulation:	0.02 M Sodium Citrate-HCl, pH 7.4		
Purity:	>95% by SDS-PAGE analysis		
Concentration:	≥ 1.0 mg/ml, lot specific		
Extinction Coefficient:	$E^{0.1\%}_{280\text{nm}} = 1.51$		
Reconstitution:	Protect from light. Thaw at 37°C without agitation until completely liquid. Then gently mix before using. Keep fibrinogen at 25-37°C. Aliquot and flash freeze unused portion.		
Storage & Stability:	Protect from light. Store at -80°C for up to 1 year. Upon initial thaw, prepare working aliquots and store at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.