

## Rabbit IgG, Protein A Purified pAb

<b>Catalog No.</b>	CSI20228A CSI20228B	<b>Quantity:</b>	500 mg 1.0 g
<b>Description:</b>	Purified from normal, non-immunized rabbits by Protein A affinity chromatography. Greater than 95% pure and preservative free.		
<b>Concentration:</b>	approximately 10 mg/ml, lot specific		
<b>Volume:</b>	Lot specific		
<b>Specificity:</b>	This antibody has been selected as an isotype control as it has no known specificity.		
<b>Source:</b>	Rabbit plasma		
<b>Molecular Weight:</b>	160 kDa		
<b>Isotype:</b>	IgG, polyclonal		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Formulation:</b>	0.01 M Sodium Phosphate, 0.15M NaCl, pH 7.4		
<b>Purification:</b>	Protein A affinity chromatography		
<b>Applications:</b>	ELISA, WB		
<b>Storage &amp; Stability:</b>	Store at -80°C product for up to 1 year. Upon initial thaw, prepare single-use aliquots and refreeze. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

