

Human Serum, Off the Clot, Pooled Donors

Catalog No.	CSI20429A	Quantity:	50 ml
	CSI20429B		500 ml

Description: Off the clot serum is produced without anticoagulants added, allowing whole blood to coagulate naturally. Serum is isolated by centrifugation.

Source: Pooled from at least 3 unique donors.

Appearance: Yellow to orange liquid

Preservatives: None

Country of Origin: USA

Storage & Stability: Store at -20°C to -80°C for at least 1 year. Upon initial thawing, prepare aliquots for storage at -20°C to -80°C. **Avoid repeated freeze/thaw cycles.**

Infectious Disease Statement: Each unit is tested for the standard FDA-required viral markers, and found negative for HBsAg, HCV, HIV-1, HIV-2, HIV-1Ag or HIV-1 NAT, ALT, West Nile Virus NAT, Zika NAT and syphilis using FDA-approved methods. While we make every effort to ensure the safety of our products, we recommend following universal precautions as laid out by the following organizations WHO, CDC and OSHA.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com