

## Mouse Anti-Human Factor X Clone 4E3F8 mAb

<b>Catalog No.</b>	CSI20515A CSI20515B	<b>Quantity:</b>	1.0 mg 10 mg
<b>Description:</b>	Mouse monoclonal antibody to human Factor X. IgG fraction purified by immobilized Protein G. Isotype IgG1.		
<b>Concentration:</b>	3.5 mg/ml		
<b>Specificity:</b>	Recognizes Factor X		
<b>Host:</b>	Mouse Full length native human purified protein.		
<b>Immunogen:</b>	Full length native protein (purified) (Human).		
<b>Isotype:</b>	IgG1		
<b>Clone:</b>	4E3F8		
<b>Formulation:</b>	Liquid in 0.05M Sodium Phosphate + 0.1M NaCl + 1 mM EDTA, pH 6.6.		
<b>Purification:</b>	IgG fraction purified by immobilized Protein G		
<b>Cross-Reactivity:</b>	Human		
<b>Applications:</b>	ELISA, Capture, WB		
<b>Storage &amp; Stability:</b>	When stored at -80°C, product is stable for 5 years from date of delivery. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.