

## US Origin Rat Plasma

|                    |           |                  |        |
|--------------------|-----------|------------------|--------|
| <b>Catalog No.</b> | CSI20577A | <b>Quantity:</b> | 1.0 ml |
|                    | CSI20577B |                  | 10 ml  |

**Description:** Plasma is prepared from healthy, non-immunized, adult Sprague Dawley rat. Blood cells are removed by centrifugation and the resulting plasma is pooled and immediately frozen.

**Specify anticoagulant:** K2EDTA, Sodium Citrate, Heparin

**Inquire about other strains:** Wistar HAN strain

**Single gender by request.**

The animals used for this material are of US Origin, ensuring that only the highest-quality material is used for the production of the plasma. We offer reliable, consistent products that deliver reliable, consistent results.

**Anticoagulant:** **Specify anticoagulant**

**Gender, Age, Strain:** Mixed gender healthy adult Sprague Dawley rat  
**Wistar HAN strain and single gender available by request.**

**Source:** Sprague Dawley rat whole blood

**Formulation:** Plasma with anticoagulant.

**Storage & Stability:** Store at -80°C as supplied for at least 1 year. Upon initial thaw, prepare working aliquots and store at -80°C. **Avoid repeated freeze-thaw cycles.**

**Country of Origin:** USA

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

