

US Origin Mouse Plasma

Catalog No.	CSI20582A	Quantity:	1.0 ml
	CSI20582B		10 ml

Description: Plasma is prepared from healthy, non-immunized, adult CD-1 mouse. Blood cells are removed by centrifugation and the resulting plasma is pooled and immediately frozen.

Specify anticoagulant: K2EDTA, Sodium Citrate, Heparin

Inquire about other strains: BALB/c and C57BL/6 strains

Inquire about single gender pooled plasma.

The animals used for this material are of US Origin, ensuring that only the highest-quality material is used for the production of the plasma. We offer reliable, consistent products that deliver reliable, consistent results.

Anticoagulant: **Specify anticoagulant.**

Gender, Age, Strain: Mixed gender healthy adult CD-1 mice
Other strains and single gender plasma available by request.

Source: Mouse whole blood

Formulation: Plasma with anticoagulant.

Storage & Stability: Store at -80°C as supplied for at least 1 year. Upon initial thaw, prepare working aliquots and store at -80°C. **Avoid repeated freeze-thaw cycles.**

Country of Origin: USA

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

