

## US Origin Mouse Serum, Pooled

<b>Catalog No.</b>	CSI20585A	<b>Quantity:</b>	1.0 ml
	CSI20585B		10 ml

**Description:** Serum is separated 'off-the-clot' from whole blood of mixed gender, adult, non-immunized, healthy CD-1 mice, pooled and immediately frozen. Sera from BALB/c, C57BL/6 and single gender pooled serum are available by request.

The animals used for this material are of US Origin, ensuring that only the highest-quality material is used for the production of the serum. We offer reliable, consistent products that deliver reliable, consistent results.

**Gender, Age, Strain:** Mixed gender, adult CD-1 mice.  
Other mouse strains and single gender serum available by request.

**Source:** US origin healthy mouse whole blood

**Formulation:** Serum, not filtered.

**Storage & Stability:** Store at -80°C as supplied for at least 1 year. Upon initial thaw, prepare working aliquots and store at -80°C. **Avoid repeated freeze-thaw cycles.**

**Country of Origin:** USA

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**