

US Origin Rat Serum, Pooled

Catalog No.	CSI20587A	Quantity:	1.0 ml
	CSI20587B		10 ml

Description: Serum is separated 'off-the-clot' from whole blood of mixed gender, adult, non-immunized, healthy Sprague Dawley (SD) rats, pooled and immediately frozen. Sera from Wistar HAN rats and single gender serum are available by request.

The animals used for this material are of US Origin, ensuring that only the highest-quality material is used for the production of the serum. We offer reliable, consistent products that deliver reliable, consistent results.

Gender, Age, Strain: Mixed gender, adult Sprague Dawley rat.
Other rat strains and single gender serum available by request.

Source: US origin healthy rat whole blood

Formulation: Serum, not filtered.

Storage & Stability: Store at -80°C as supplied for at least 1 year. Upon initial thaw, prepare working aliquots and store at -80°C. **Avoid repeated freeze-thaw cycles.**

Country of Origin: USA

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.