

Research Grade Pagibaximab

Catalog No.	DGX03201A DGX03201B	Quantity:	100 µg 1.0 mg
Alternate Names:	BSYX-A110, PMAB, CAS: 595566-61-3		
Description:	Research Grade Pagibaximab is a chimeric monoclonal antibody directed against <i>Staphylococcus epidermidis</i> Lipoteichoic Acid		
Target:	<i>Staphylococcus epidermidis</i> Lipoteichoic Acid		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Chimeric IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Pagibaximab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

