

TNFRSF10B Research Grade Tigatuzumab

Catalog No.	DHA29102A DHA29102B	Quantity:	100 µg 1.0 mg
Alternate Names:	CS-1008, TRA-8, CAS: 918127-53-4		
Description:	Research Grade Tigatuzumab is a humanized monoclonal antibody which targets human Tumor necrosis factor receptor superfamily member 10B (TRAIL-R2, CD262)		
UniProt ID (target):	O14763		
Specificity:	Human TNFRSF10bB (CD262)		
Source:	XtenCHO cell line		
Isotype:	Human IgG1-nd		
INN:	Tigatuzumab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.