

TNFSF11 Research Grade Denosumab

Catalog No.	DHA30301A DHA30301B	Quantity:	100 µg 1.0 mg
Alternate Names:	AMG162, CAS: 615258-40-7		
Description:	Research Grade Denosumab is a human monoclonal antibody directed against Human Tumor necrosis factor ligand superfamily member 11 (RANKL, CD254)		
UniProt ID (target):	O14788		
Target:	Human TNFSF11 (RANKL, CD254)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG2 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Denosumab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.