

TNFSF14

Research Grade Quisovalimab

Catalog No.	DHA70801A DHA70801B	Quantity:	100 µg 1.0 mg
Alternate Names:	AEVI-002, CERC-002, F-19, MDGN-002, SAR-252067 CAS: 2427667-03-4		
Description:	Research Grade Quisovalimab is a human monoclonal antibody that targets human Tumor necrosis factor superfamily member 14, Herpes virus entry mediator ligand (HVEM-L) CD258		
UniProt ID (target):	O43557		
Specificity:	Human TNFSF14 (HVEM-L, CD258)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO cell line		
Isotype:	Human IgG4 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
Applications:	Functional studies		
INN:	Quisovalimab		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

