

PSCA

Research Grade MK-4721

Catalog No.	DHA73501A DHA73501B	Quantity:	100 µg 1.0 mg
Alternate Names:	AGS-PSCA, MK-4721		
Description:	Research Grade MK-4721 is a human monoclonal antibody which targets Human Prostate stem cell antigen		
UniProt ID (target):	O43653		
Specificity:	Human Prostate stem cell antigen		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A or G affinity chromatography		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

