

TNFRSF13B-Fc Research Grade Atacicept

Catalog No.	DHB34103A DHB34103B	Quantity:	100 µg 1.0 mg
Alternate Names:	TNFRSF13B-Fc, TACI-Fc, CAS: 845264-92-8		
Description:	Atacicept is a fusion protein combining TNFRSF13B and Fc (fragment constant) of IgG1 heavy chain		
UniProt ID:	O75888 & Q9Y275		
Source:	Mammalian cell line		
Concentration:	1.0 mg/ml		
Endotoxin:	< 0.01 EU/µg by LAL analysis		
Formulation:	0.01M PBS, pH 7.4		
Purity:	> 95% as determined by SDS-PAGE		
Storage & Stability:	Store at 2-8°C for 1 week or as aliquots for up to 1 year at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.