

DKK1

Research Grade Anti-Human DKK1 (BHQ880) mAb

Catalog No.	DHB48102A DHB48102B	Quantity:	100 µg 1.0 mg
Alternate Names:	BHQ880, BHQ-880, BHQ 880, Dickkopf-related protein 1, DKK1, SK, Dickkopf-1, Dkk-1, hDkk-1		
Description:	Research Grade Anti-Human DKK1 (BHQ880) is human monoclonal antibody that binds to and inhibits human Dickkopf-related protein 1 (Dkk-1)		
UniProt ID (target):	O94907		
Specificity:	Human Dickkopf-related protein 1 (Dkk-1)		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Human IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	Functional studies, Neutralizing		
Clone:	BHQ880		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

