

LYPD3

Research Grade Lupartumab

Catalog No.	DHB60001A DHB60001B	Quantity:	100 µg 1.0 mg
Alternate Names:	BAY 1112623, BAY 1129980, CAS: 1640971-88-5		
Description:	Research Grade Lupartumab is a human monoclonal antibody which targets human Ly6/PLAUR domain-containing protein 3 (MIG-C4)		
UniProt ID (target):	O95274		
Target:	Human Ly6/PLAUR domain-containing protein 3		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG1 lambda		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Lupartumab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

