

NTN1

Research Grade Anti-Human NTN1/Netrin-1 (NP137) mAb

Catalog No.	DHB70501A DHB70501B	Quantity:	100 µg 1.0 mg
Alternate Names:	Netrin-1, NTN1, Epididymis tissue protein Li 131P, NTN1L		
Description:	Research Grade Anti-Human NTN1/Netrin-1 (NP137) is a human monoclonal which targets human Netrin-1.		
UniProt ID (target):	O95631		
Specificity:	Human Netrin-1		
Source:	XtenCHO cell line		
Isotype:	Human IgG1-kappa		
Clone:	NP137		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4.		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com