

IL33

Research Grade Torudokimab

Catalog No.	DHB74003A DHB74003B	Quantity:	100 µg 1.0 mg
Alternate Names:	LY3375880, CAS: 2241728-76-5		
Description:	Research Grade Torudokimab is a humanized monoclonal antibody directed against human Interleukin-33		
UniProt ID (target):	O95760		
Target:	Human Interleukin-33		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG4 kappa		
Formulation:	Sterile-filtered 0.01 M PBS, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Torudokimab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.