

LY6G6D

Research Grade Anti-Human LY6G6D/MEGT1 (1G4) mAb

Catalog No.	DHB78001A DHB78001B	Quantity:	100 µg 1.0 mg
Alternate Names:	Lymphocyte antigen 6 complex locus protein G6d, Protein Ly6-D, Megakaryocyte-enhanced gene transcript 1 protein		
Description:	Research Grade Anti-Human LY6G6D/MEGT1 (1G4) mAb is a humanized monoclonal antibody which targets human Protein Ly6-D		
UniProt ID (target):	O95868		
Specificity:	Human Protein Ly6-D		
Source:	Mammalian cells		
Isotype:	Human IgG1		
Clone:	1G4		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4.		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Application:	Functional studies		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com